

## CERTIFICATE OF TYPE TESTING IN ACCORDANCE WITH EN 14175

CERTIFICATE & REPORT NO: INV/EN14175/568-A  
ISSUE DATE: 24th April 2013

**Fume Cupboard Type:**

fumetop1800

**FH-180-PP bench type**

*External Dimensions:*

Height = 2370mm  
(1500mm + 870mm)

Width = 1800mm

Depth = 850mm

*Internal Dimensions:*

Height = 1190mm

Width = 1790mm

Depth = 620mm (wall-sash)  
575mm (baffle-sash)

*Baffle Gap Dimensions:*

Top gap: 40mm

Bottom gap: 230mm

**Fume Cupboard Manufacturer:**

**TopAir USA**

**3182 Monterey Drive**

**Merrick**

**NY 11566**

**USA**



**Fume Cupboard Flow:**

**Volume flow rate: 0.470m<sup>3</sup>/s**

**Test opening width: 1610mm**

**Test opening height: 500mm**

**Face velocity: 0.515m/s**

**Pressure drop: 125Pa @ 250mmD**

This is to certify that the fume cupboard described above has been type tested in accordance with Part 3, in compliance with the requirements of Part 2 and with reference to Part 1 of EN 14175 and resulted in performance characteristics given in the corresponding test report.

**Tested and Certified by:** Dr A F Bicen

Invent UK Ltd, 85 Southdown Road, Harpenden, Herts AL5 1PR

Tel: 01582 461666 : Fax: 01582 460077 : E-mail: info@invent-uk.com

A handwritten signature in blue ink that reads 'A.F. Bicen'.