

CERTIFICATE OF TYPE TESTING IN ACCORDANCE WITH EN 14175

CERTIFICATE & REPORT NO: INV/EN14175/558-A
ISSUE DATE: 24th April 2013

Fume Cupboard Type:

fumetop1800

FH-180 bench type

External Dimensions:

Height = 2300mm
(1500mm + 800mm)

Width = 1800mm

Depth = 800mm

Internal Dimensions:

Height = 1050mm

Width = 1550mm

Depth = 670mm (wall-sash)
620mm (baffle-sash)

Baffle Gap Dimensions:

Top gap: 40mm

Bottom gap: 140mm - 240mm

Fume Cupboard Manufacturer:

TopAir USA

3182 Monterey Drive

Merrick

NY 11566

USA



Fume Cupboard Flow:

Volume flow rate: 0.430m³/s

Test opening width: 1550mm

Test opening height: 500mm

Face velocity: 0.516m/s

Pressure drop: 95Pa @ 250mmD

This is to certify that the fume cupboard described above has been type tested in accordance with Part 3, in compliance with the requirements of Part 2 and with reference to Part 1 of EN 14175 and resulted in performance characteristics given in the corresponding test report.

Tested and Certified by: Dr A F Bicen

Invent UK Ltd, 85 Southdown Road, Harpenden, Herts AL5 1PR

Tel: 01582 461666 : Fax: 01582 460077 : E-mail: info@invent-uk.com

A handwritten signature in blue ink that reads 'A.F. Bicen'.